

SURGICAL DRESSINGS

ISLAND DRESSINGS - WATERPROOF



Application:



1. Peel away both paper backing sheets on reverse side of the dressing.
2. Apply dressing to the wound site and remove large fixation sheet, slowly, parallel to the skin in the direction of arrows shown.
3. Remove remaining edging paper strip in direction of arrows.

CODE	ADHESIVE SHEET SIZE	ABSORBENT PAD SIZE	PACKAGING
SED006	38 x 38mm	20 x 20mm	Box 100
SED007	40 x 70mm	20 x 40mm	Box 100
SED008	60 x 60mm	30 x 30mm	Box 100
SED009	60 x 90mm	30 x 60mm	Box 50
SED010	90 x 150mm	50 x 100mm	Box 30
SED011	110 x 210mm	50 x 150mm	Box 20

PRODUCT NAME

Sere-derm Plus

DESCRIPTION

Sterile transparent waterproof dressing with non-adhesive absorbent pad

FEATURES

- Small sized **Bosui Pad**
- Reduces risk of secondary infection.
- Protects wound from penetration of bacteria.
- Waterproof - allows bathing etc.
- Allows skin to breath, prevents maceration.
- Adequate strength preventing breakage of film when removed.

USES

- Ideal for minor cuts and wounds
- Minor operative procedure sites - i.e. face.
- Following skin biopsy sites.
- Hormonal implant and antitubercle sites.
- For use after skin lesion treatment with **CryoPen**.

Sere-med Dressings - *The Affordable Choice*