

**SURGICAL DRESSINGS**

**ROLL DRESSING - FABRIC**



**PRODUCT NAME**

**Sere-Fix Roll**

**DESCRIPTION**


Non-Sterile fabric. Latex-free adhesive surgical dressing on a roll.

**FEATURES**

- Soft non-woven fabric.
- Latex free.
- Hypoallergenic adhesive.
- Cut to size required affording smaller or larger wound sites to be dressed\*.
- No need to stock a variety of surgical dressings of different sizes.
- Can be used on a variety of different sized wounds / injuries\*.
- Affords accessibility and dressing of difficult wound sites\*.

**USES**

- Dressing of wounds.
- Also used as a protective layer on the skin prior to the use of strapping or casting products etc.



\* Use with Sere-pad instead of traditional gauze to create a surgical dressing with an absorbent pad.

CODE	SIZE	PACKAGING
SER 017	100mm x 10m	Box 10m

**Sere-med Dressings - *The Affordable Choice***