

**SURGICAL DRESSINGS**

**NON-STICK DRESSING**



Please note dressing supplied is a white absorbent sheet.

**PRODUCT NAME**

Sere-Pad

**DESCRIPTION**

Sterile non-stick dressing is a highly absorbable pad protected on both sides by a thin but durable film covering.

**FEATURES**

- Pad does not stick to the wound.
- Cut to size required affording smaller wound sites to be dressed.
- Absorbs wound exudate.
- Affords accessibility and dressing of difficult wound sites\*.
- Keeps wound area dry.

**USES**

- Use instead of traditional gauze when dressing, for example, exudating wounds.

\* Use Sere-pad instead of traditional gauze with SurginLeep (sterile) / Sere-derm Roll Film Dressings or Sere-Fix Fabric Roll to create either a waterproof or a surgical dressing with an absorbent pad.

*Compare to Telfa Dressing*

CODE	SIZE	PACKAGING
SER 011	60 x 100mm	Box 50

**Sere-med Dressings - *The Affordable Choice***