


SURGICAL DRESSINGS

ISLAND DRESSINGS - FABRIC



PRODUCT NAME

SurginPore

DESCRIPTION

Sterile surgical adhesive plaster dressing with a non-adhesive absorbable pad

FEATURES

- Non-woven fabric.
- Innovative application of adhesive assures superior air permeability affording the skin to breathe.
- Water repellent.
- Soft fabric - ideal for use on elderly patients or sensitive skins.
- Hypo-allergenic adhesive.

USES

- Post operative care where 12/24 hourly dressings are required.
- Minor procedures where long term (**waterproof**) dressings are not required.

Compare to Primapore

Pharmacy Retail Display Box
Available

CODE	ADHESIVE SHEET SIZE	ABSORBENT PAD SIZE	PACKAGING
SUR001	50 x 80mm	20 x 40mm	Box 50
SUR002	60 x 100mm	30 x 60mm	Box 50
SUR003	80 x 130mm	40 x 90mm	Box 30
SUR004	100 x 160mm	50 x 100mm	Box 30
SUR006	110 x 200mm	50 x 150mm	Box 20
SUR007	120 x 260mm	60 x 200mm	Box 20