

# Effective Application Times for The Veterinary Practice



CryoPen®



**Sere-med**  
Life Enhancing Innovations



### Perianal Adenoma

Application 15-20 seconds 2 applications per sitting, follow up treatment required.



### Carcinoma

Application 10 seconds 2 applications per sitting, follow up treatment required

<https://www.ho-equipments.com/practitioners/veterinary/>

**For video material please use above link.**



### Cyst

Application 10 seconds 2 applications per sitting, follow up treatment required.



### Papilloma

Application 10-15 seconds 2 applications per sitting, follow up treatment required.

Application 6 seconds 1 application per sitting.

### Skin Tags



**Sere-med**  
Life Enhancing Innovations