

LUXAMED®

Medizintechnik Made in Germany
Specializing in Accuracy

Precision for the Professional

Introducing the:
Blood Pressure Meter Precision Dial PD

The LUXAMED® Dial type blood pressure meter guarantees quality, reliability, highest precision, user friendliness and provides a unique colour coded measurement range category.

HOW THIS FEATURE WORKS

During blood pressure measurement, the colour on the dial will indicate the range category in which the reading was observed. This assists the Practitioner to easily take an accurate reading and diagnose the extent of possible hypertension.

Measurement Category Range Table:

Colour	Category	Systolic Range	Diastolic Range	Description
Green	Normotensive	≤ 120	≤ 80	Normal BP
Yellow	Prehypertension	120 - 139	80 - 90	At risk / Elevated BP
Red	Hypertension	>140	>90	High BP

MORE ABOUT THE PD

The LUXAMED® precision dial PD Blood Pressure Meter has a non-corrosive precision measurement mechanism encased within a shock protective housing. This aneroid dial type BP meter utilizes a two tube and inflation bulb connection system that allows for a **one hand operation**. A 60mm above standard sized dial face ensures the measurement information is clear, easily visible and can be accurately observed.

Supplied in a portable zip bag, the precision dial is supplied with a 22-32 cm standard adult inflation system.

OPTIONAL ACCESSORIES

*Child Inflation System (Velcro cuff and inflation bag)
Size 14cm to 21cm*



*Large Adult Inflation System (Velcro cuff and inflation bag)
Size 33cm to 41cm*



Orders ☎ +27 11 262 2220



www.seremed.com

SAHPRA Licensed