



Sere-med Diagnostic LED Penlight



LUXAMED® DIAGNOSTICS Clarity in the right light

High-quality, high-power LED Penlight for eye and throat diagnostic.

Characteristics of our LUXAMED® Diagnostic LED Penlight

- Available in red or blue.
- Best quality ABS plastic housing
- Luminous intensity: 20 000 – 25 000 mcd (mcd is the standard unit for measuring the intensity of a light source).
- Lifetime of the LED – 100 000 hours
- Batteries (2 x AAA) – 50 hours illumination time.
- Operating voltage – 3V
- Weight – 31.5g
- Fingertip on / off clip type release switch

Uses:

Throat

To identify an infection or inflammation using white light.

Eyes

Important – LED is approved for eye diagnostic

- Pupil dilation

Seeing the difference with the Blue Cobalt filter

Place the blue cobalt filter over the front end of the LED bulb. By using fluorescence in the eye, the blue light will illuminate scratching and abrasions more clearly.



Orders ☎ +27 11 262 2220



www.seremed.com

SAHPRA Licensed

LUXAMED®
Medizintechnik Made in Germany