



Sere-med OTOSCOPE



An Otoscope affording an Unrestricted Field of View from Luxamed,[®] Germany

CCT - Colour changing Technology

A maintenance-free lightweight pocket sized 2.5V LED CCT - Colour changing Technology Otoscope offering the following features:

1. One switch operation with automatic switch-off after 3 minutes saves battery life.
2. Optimized head design with insufflation port (sufflation device is optional and not included) and a lockable pivoting magnifying glass, that when moved to the left side, affords the use of the appropriate ear instrumentation.
3. An aluminium fibre reinforced handle for 2 x AAA alkaline batteries.
4. No bulb replacement, therefore, no additional costs as Luxamed Otoscope LEDs have a life expectancy of 100 000 hours due to a constant current drive which does not change the intensity overtime, therefore prolonging the life of the LED bulbs.
5. Supplied with 20 grey ear speculae (10 x size 2.5 mm) and (10 x size 4.0 mm).

Why is the Luxamed CCT Otoscope Unique?

Superior illumination is guaranteed by LED bulbs that are concentrically placed within the Otoscope head.

This is a key feature of the Luxamed Otoscope where light is shone from every possible angle thereby eliminating any outside factors of light that may influence what the medical practitioner is looking at.

CCT 2.5 LED

- Innovative illumination provided by 2 x 4 LEDs in a concentric ring design.
- Depressing the on / off switch which affords light options of either:



WARM YELLOW LIGHT

Warm light (colour temperature 3000 k / warm yellowish colour) which helps to see the bone structure in the ear better, liquid behind the tympanic membrane (eardrum) or identifying a ruptured or torn eardrum with a little more certainty.



COLD WHITE LIGHT

Cold Light (colour temperature 6000 k / cold white blueish colour) ideally suited for looking at ear inflammation and infection in the beginning stages of these disorders. Under White light, areas with high blood flow (capillaries) make a better contrast to the white light which illuminates ear infection and inflammation more clearly. Infections under white light can be detected during the earlier stages.

Orders  +27 11 262 2220



www.seremed.com

SAHPRA Licensed

LUXAMED[®]
Medizintechnik Made in Germany