



Sere-med STETHOSCOPE



LUXAMED® Germany, introduces a game changer in diagnostic acoustics

STETHOSCOPE CX CARDIOLOGY

"I have never heard such QUALITY from a stethoscope"

LUXAMED® Cardiology stethoscopes are characterized by a high mass CNC-milled manufactured processed stainless steel chest piece. This produces a high inertia which reduces vibrations of the chest piece resulting in significantly improved audibility in all frequency ranges.

STETHOSCOPE CARDIOLOGY CX therefore offers excellent sound characteristics due to the perfectly tuned mass proportion of the reversible chest piece to the vibration diaphragms.

BINAURAL

The binaural of the Cardiology CX is 72cm in length and longer than most Cardiology stethoscopes. The additional length affords a more comfortable working distance without compromising sound quality.

THE BELL

The standard bell configuration with diaphragm, is ideal for use on smaller parts of the body. To enhance the function of the bell, simply remove the small vibration diaphragm to listen to lower frequency sounds.

THE RETAINING RINGS AND DIAPHRAGMS

These are made of light metal with a SOFT FEEL coating for more comfort on the skin.

The retaining rings and diaphragms are easy to disinfect, guarantee limited wear and tear and the chance of breaking is reduced. The retaining rings facilitates easy exchanging of both vibration diaphragms.

The chest piece (diaphragm and bell sections) have very fine threads which reduces dirt collection and ensure minimal cleaning. The design facilitates easy removal, replacement and or re-attachment of the retaining rings and diaphragms especially following the use of the bell in the low frequency application.

SOFT SILICONE REMOVABLE EAR TIPS

Adapt optimally to the auditory canal affording a comfortable fit in the ear.

These unique SOFTBUD Ear Tips in a sound insulating design, guarantee optimum sound transmission and perfect auscultation without any background noise. The inside plastic body of the Ear Tips prevent the inadvertent closing of the Ear Tips during head movements and ensures a noiseless movement of the binaural.

SOFTBUD Ear Tips – small, medium, and large supplied with the CX are easy to clean, remove, re-attach and replace.

SPECIFICATIONS

Latex and nickel free binaural, chest piece and Ear Tips.

Weight	(Chest piece): 160gm
Diameter	(Large Diaphragm): 3.9cm
Diameter	(Small Diaphragm): 3.2cm
Length	(Binaural incl Ear tips): 72cm
Length	(Total incl Chest piece): 77cm

LUXAMED®
Medizintechnik Made in Germany

Orders  +27 11 262 2220



www.seremed.com

SAHPRA Licensed